PROSTAT® CORPORATION



CPM-760 6inch Charge Plate Assembly

The Prostat CPM-760A is a 20pF, 6-inch Charge Plate Monitor assembly intended for use with the Prostat PFM-711A Field Meter, PCS-730 Charging Source and PDT-740B Static Decay Timer. This handy accessory is designed to perform very accurate performance audits on overhead ionizers and room ionization systems. Used with extension wings, the CPM-760A is also used to evaluate and balance bench top ionizers.

The CPM-760A is completely collapsible for compact storage in the Prostat Professional ESD Auditing Kit Series. It may be used for evaluating voltage generation and static decay performance of personnel, footwear, flooring, chairs, carts and other equipment and materials.

Specifications

| Isolated Plate: | Electro-plated Aluminum plate, 6.0" square |
|-----------------|--|
| Ground Plate: | Plated aluminum sheet, 6.0" square |
| Plate Spacing: | Approximately 0.63" between Isolated and Ground plates. Plate spacers made from machined virgin, white Teflon® Rod |
| Support: | Stored disassembled, plated aluminum sheet and steel lower stabilizing block |
| Assembly: | Assembled with machine cap screws |
| Capacitance: | Parallel Plate Capacitance 20 picofarads (pF) ±2 pF |
| Range: | 0 to \pm > 5 kV at a plate separation of 0.63" |
| Accuracy: | ±2% 0 to 1,000 Volts |

| | 5 Minute Decay when charged to 1kV, floating plate voltage decays less than 5% (< 50 @ 50% Rh: volts) over a 5 minute period |
|------------------|--|
| Calibration: | Accomplished via adjustable instrument saddle and grounding device |
| Grounding: | Via snap fastener of Model PFM-711A and meter saddle assembly |
| Assembly Height: | 6.0" High in Vertical Test Configuration 9.0" High in Horizontal Test Configuration; Center of Plate 6.0" from Worksurface |
| Weight: | 2.5 lbs (1.13 kg) |



©2007 Prostat Corporation. All rights reserved. Prostat is the registered trademarks of Prostat Corporation in the United States and other countries.

1072 Tower Lane, Bensenville, IL 60106 USA Phone +1 (630) 238-8883 Fax +1 (630) 238-9717 Web access: http://www.prostatcorp.com

Specifications subject to change without notice. Printed in U.S.A. REV2: 6/20/07



Pro-Pack Materials Pte. Ltd. Blk 970, Toa-Payoh North #04-18/21 Toa Payoh Industrial Estate, Singapore 318992 (65) 6-354-1511 Fax: (65) 6-354-1711

E-mail: sales@pro-pack.com.sg Web access: http://www.pro-pack.com.sg